

Product datasheet for **GA117403**

Human TMEM95 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM95 CRISPRa kit – CRISPR gene activation of human transmembrane protein 95
Symbol:	TMEM95
Components:	GA117403G1, TMEM95 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117403G2, TMEM95 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117403G3, TMEM95 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q3KNT9
Synonyms:	UNQ9390
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

