

Product datasheet for **GA117421**

Human LINGO4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LINGO4 CRISPRa kit - CRISPR gene activation of human leucine rich repeat and Ig domain containing 4
Symbol:	LINGO4
Components:	GA117421G1, LINGO4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117421G2, LINGO4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117421G3, LINGO4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6UY18</u>
Synonyms:	DAAT9248; LRRN6D; PRO34002
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)