

## Product datasheet for **GA117428**

### Human FAM5C (BRINP3) activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | BRINP3 CRISPRa kit - CRISPR gene activation of human BMP/retinoic acid inducible neural specific 3 |
| Symbol:       | BRINP3   |
| UniProt ID:   | <a href="#">Q76B58</a>   |
| Synonyms:     | DBCCR1L; DBCCR1L1; FAM5C   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                            |

#### Product images:



[View online >](#)