

Product datasheet for **GA117430**

Human TFAP2E activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TFAP2E CRISPRa kit - CRISPR gene activation of human transcription factor AP-2 epsilon
Symbol:	TFAP2E
Components:	GA117430G1, TFAP2E gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117430G2, TFAP2E gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117430G3, TFAP2E gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6VUC0</u>
Synonyms:	AP2E
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)