

Product datasheet for **GA117432**

Human Marapsin2 (PRSS38) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PRSS38 CRISPRa kit - CRISPR gene activation of human serine protease 38
Symbol:	PRSS38
Components:	GA117432G1, Marapsin2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117432G2, Marapsin2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117432G3, Marapsin2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A1L453</u>
Synonyms:	MPN2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)