


Product datasheet for **GA117437**

Human EAN57 (TEX33) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TEX33 CRISPRa kit - CRISPR gene activation of human testis expressed 33
Symbol:	TEX33
Components:	GA117437G1, EAN57 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117437G2, EAN57 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117437G3, EAN57 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>O43247</u>
Synonyms:	C22orf33; cE81G9.2; EAN57
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)