

Product datasheet for **GA117443**

Human C2orf53 (PRR30) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PRR30 CRISPRa kit - CRISPR gene activation of human proline rich 30
Symbol:	PRR30
Components:	GA117443G1, C2orf53 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117443G2, C2orf53 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117443G3, C2orf53 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q53SZ7</u>
Synonyms:	C2orf53
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

