

Product datasheet for **GA117444**

Human C2orf74 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	C2orf74 CRISPRa kit - CRISPR gene activation of human chromosome 2 open reading frame 74
Symbol:	C2orf74
Components:	GA117444G1, C2orf74 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117444G2, C2orf74 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117444G3, C2orf74 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A8MZ97</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

