


## Product datasheet for **GA117450**

### Human **CCDC158** activation kit by **CRISPRa**

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                  |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                     |
| Description:  | CCDC158 CRISPRa kit - CRISPR gene activation of human coiled-coil domain containing 158                                                                                                                                                                                                                          |
| Symbol:       | CCDC158                                                                                                                                                                                                                                                                                                          |
| Components:   | GA117450G1, CCDC158 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA117450G2, CCDC158 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA117450G3, CCDC158 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q5M9N0</a></u>                                                                                                                                                                                                                                                                                    |
| Synonyms:     | FLJ25770; MGC35086                                                                                                                                                                                                                                                                                               |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                          |

#### Product images:



[View online »](#)