

## Product datasheet for **GA117458**

### Human ARSI activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	ARSI CRISPRa kit - CRISPR gene activation of human arylsulfatase family member I
Symbol:	ARSI
UniProt ID:	<a href="#">Q5FYB1</a>
Synonyms:	ASI; SPG66
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)