

Product datasheet for GA117460

Human SLC35D3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SLC35D3 CRISPRa kit - CRISPR gene activation of human solute carrier family 35 member D3
Symbol:	SLC35D3
UniProt ID:	<u>Q5M8T2</u>
Synonyms:	FRCL1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)