

Product datasheet for **GA117465**

Human TREML1 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | TREML1 CRISPRa kit - CRISPR gene activation of human triggering receptor expressed on myeloid cells like 1 |
| Symbol: | TREML1 |
| UniProt ID: | Q86YW5 |
| Synonyms: | dj238O23.3; GLTL1825; PRO3438; TLT-1; TLT1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)