

Product datasheet for **GA117472**

Human Tetraspanin 33 (TSPAN33) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TSPAN33 CRISPRa kit - CRISPR gene activation of human tetraspanin 33
Symbol:	TSPAN33
UniProt ID:	Q86UF1
Synonyms:	PEN; PEN.; TSPAN-33
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)