

Product datasheet for **GA117478**

Human TMEM249 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM249 CRISPRa kit - CRISPR gene activation of human transmembrane protein 249
Symbol:	TMEM249
Components:	GA117478G1, TMEM249 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117478G2, TMEM249 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117478G3, TMEM249 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q2WGI8</u>
Synonyms:	C8ORFK29
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)