

Product datasheet for **GA117493**

Human RTL4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	RTL4 CRISPRa kit - CRISPR gene activation of human retrotransposon Gag like 4
Symbol:	RTL4
Components:	GA117493G1, RTL4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117493G2, RTL4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117493G3, RTL4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6ZR62</u>
Synonyms:	Mar4; Mart4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)