

Product datasheet for **GA117508**

Human MS4A10 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MS4A10 CRISPRa kit - CRISPR gene activation of human membrane spanning 4-domains A10
Symbol:	MS4A10
UniProt ID:	Q96PG2
Synonyms:	CD20L7; MS4A9
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)