

Product datasheet for **GA117517**

Human Ankyrin repeat domain containing protein 33 (ANKRD33) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ANKRD33 CRISPRa kit - CRISPR gene activation of human ankyrin repeat domain 33
Symbol:	ANKRD33
Components:	GA117517G1, Ankyrin repeat domain containing protein 33 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117517G2, Ankyrin repeat domain containing protein 33 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117517G3, Ankyrin repeat domain containing protein 33 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q7Z3H0
Synonyms:	C12orf7; PANKY
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)