

Product datasheet for **GA117525**

Human SLC35F4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SLC35F4 CRISPRa kit – CRISPR gene activation of human solute carrier family 35 member F4
Symbol:	SLC35F4
Components:	GA117525G1, SLC35F4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117525G2, SLC35F4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117525G3, SLC35F4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A4IF30</u>
Synonyms:	C14orf36; c14_5373; FLJ37712
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

