

## Product datasheet for **GA117530**

### Human TMC3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TMC3 CRISPRa kit - CRISPR gene activation of human transmembrane channel like 3
Symbol:	TMC3
Components:	GA117530G1, TMC3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117530G2, TMC3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117530G3, TMC3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u><a href="#">Q7Z5M5</a></u>
Synonyms:	transmembrane channel-like 3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)