

Product datasheet for **GA117534**

Human ZKSCAN2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZKSCAN2 CRISPRa kit - CRISPR gene activation of human zinc finger with KRAB and SCAN domains 2
Symbol:	ZKSCAN2
Components:	GA117534G1, ZKSCAN2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117534G2, ZKSCAN2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117534G3, ZKSCAN2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q63HK3</u>
Synonyms:	ZNF694; ZSCAN31; ZSCAN34
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

