

Product datasheet for **GA117536**

Human PKDIL3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PKDIL3 CRISPRa kit – CRISPR gene activation of human polycystin 1 like 3, transient receptor potential channel interacting
Symbol:	PKDIL3
Components:	GA117536G1, PKDIL3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117536G2, PKDIL3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117536G3, PKDIL3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q7Z443</u>
Synonyms:	polycystic kidney disease 1 like 3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

