

Product datasheet for **GA117537**

Human CD300E activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CD300E CRISPRa kit - CRISPR gene activation of human CD300e molecule
Symbol:	CD300E
UniProt ID:	Q496F6
Synonyms:	CD300LE; CLM-2; CLM2; CMRF35-A5; IREM-2; IREM2; PlgR-2; PlgR2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)