


Product datasheet for **GA117547**

Human Syncollin (SYCN) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SYCN CRISPRa kit - CRISPR gene activation of human syncollin
Symbol:	SYCN
Components:	GA117547G1, Syncollin gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117547G2, Syncollin gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117547G3, Syncollin gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q0VAF6</u>
Synonyms:	INSSA1; SYL
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)