

Product datasheet for **GA117549**

Human ZNF404 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF404 CRISPRa kit – CRISPR gene activation of human zinc finger protein 404
Symbol:	ZNF404
Components:	GA117549G1, ZNF404 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117549G2, ZNF404 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117549G3, ZNF404 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q494X3
Synonyms:	MGC120455
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

