

Product datasheet for **GA117573**

Human KCNT2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KCNT2 CRISPRa kit - CRISPR gene activation of human potassium sodium-activated channel subfamily T member 2
Symbol:	KCNT2
UniProt ID:	<u>Q6UVM3</u>
Synonyms:	KCa4.2; SLICK; SLO2.1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)