

Product datasheet for **GA117574**

Human **BARHL2** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	BARHL2 CRISPRa kit - CRISPR gene activation of human BarH like homeobox 2
Symbol:	BARHL2
Components:	GA117574G1, BARHL2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117574G2, BARHL2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117574G3, BARHL2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NY43</u>
Synonyms:	BarH-like 2; BarH-like homeobox 2; OTTHUMP00000011704
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)