

Product datasheet for **GA117590**

Human SH3RF3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SH3RF3 CRISPRa kit - CRISPR gene activation of human SH3 domain containing ring finger 3
Symbol:	SH3RF3
UniProt ID:	Q8TEJ3
Synonyms:	POSH2; SH3MD4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)