

Product datasheet for **GA117597**

Human TMPRSS7 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMPRSS7 CRISPRa kit - CRISPR gene activation of human transmembrane serine protease 7
Symbol:	TMPRSS7
Components:	GA117597G1, TMPRSS7 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117597G2, TMPRSS7 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117597G3, TMPRSS7 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q7RTY8</u>
Synonyms:	B230219I23Rik; Gm1748; matriptase-3; transmembrane serine protease 7
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)