

Product datasheet for **GA117598**

Human CD200R1L activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CD200R1L CRISPRa kit - CRISPR gene activation of human CD200 receptor 1 like
Symbol:	CD200R1L
Components:	GA117598G1, CD200R1L gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117598G2, CD200R1L gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117598G3, CD200R1L gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6Q8B3</u>
Synonyms:	CD200R2; CD200RLa
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)