

Product datasheet for **GA117599**

Human PAQR9 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PAQR9 CRISPRa kit - CRISPR gene activation of human progesterin and adipoQ receptor family member 9
Symbol:	PAQR9
Components:	GA117599G1, PAQR9 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117599G2, PAQR9 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117599G3, PAQR9 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6ZVX9</u>
Synonyms:	FLJ41938
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)