

Product datasheet for **GA117606**

Human LRIT3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LRIT3 CRISPRa kit - CRISPR gene activation of human leucine rich repeat, Ig-like and transmembrane domains 3
Symbol:	LRIT3
Components:	GA117606G1, LRIT3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117606G2, LRIT3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117606G3, LRIT3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3SXY7</u>
Synonyms:	CSNB1F; FIGLER4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)