

## Product datasheet for **GA117610**

### Human PFN3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PFN3 CRISPRa kit - CRISPR gene activation of human profilin 3
Symbol:	PFN3
UniProt ID:	<a href="#">P60673</a>
Synonyms:	1700012P12Rik; profilin 3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)