

Product datasheet for **GA117612**

Human PLCXD3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLCXD3 CRISPRa kit - CRISPR gene activation of human phosphatidylinositol specific phospholipase C X domain containing 3
Symbol:	PLCXD3
Components:	GA117612G1, PLCXD3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117612G2, PLCXD3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117612G3, PLCXD3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q63HM9
Synonyms:	phosphatidylinositol-specific phospholipase C, X domain containing 3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)