

Product datasheet for **GA117621**

Human EYS activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	EYS CRISPRa kit - CRISPR gene activation of human eyes shut homolog
Symbol:	EYS
Components:	GA117621G1, EYS gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117621G2, EYS gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117621G3, EYS gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5T1H1</u>
Synonyms:	bA74E24.1; bA166P24.2; bA307F22.3; C6orf178; C6orf179; C6orf180; dj22i17.2; dj1018A4.2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)