

Product datasheet for **GA117630**

Human **MOGAT3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	MOGAT3 CRISPRa kit - CRISPR gene activation of human monoacylglycerol O-acyltransferase 3
Symbol:	MOGAT3
Components:	GA117630G1, MOGAT3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117630G2, MOGAT3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117630G3, MOGAT3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q86VF5</u>
Synonyms:	DC7; DGAT2L2; MGAT3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)