

Product datasheet for **GA117654**

Human ARSH activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ARSH CRISPRa kit - CRISPR gene activation of human arylsulfatase family member H
Symbol:	ARSH
Components:	GA117654G1, ARSH gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117654G2, ARSH gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117654G3, ARSH gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5FYA8</u>
Synonyms:	sulfatase
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)