

Product datasheet for **GA117659**

Human LRRTM3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LRRTM3 CRISPRa kit - CRISPR gene activation of human leucine rich repeat transmembrane neuronal 3
Symbol:	LRRTM3
Components:	GA117659G1, LRRTM3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117659G2, LRRTM3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117659G3, LRRTM3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q86VH5</u>
Synonyms:	MGC131810
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)