

Product datasheet for **GA117663**

Human SERINC2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SERINC2 CRISPRa kit - CRISPR gene activation of human serine incorporator 2
Symbol:	SERINC2
Components:	GA117663G1, SERINC2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117663G2, SERINC2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117663G3, SERINC2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96SA4</u>
Synonyms:	FKSG84; PRO0899; TDE2; TDE2L
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)