

Product datasheet for **GA117679**

Human SELV (SELENOV) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SELENOV CRISPRa kit - CRISPR gene activation of human selenoprotein V
Symbol:	SELENOV
Components:	GA117679G1, SELV gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117679G2, SELV gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117679G3, SELV gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P59797</u>
Synonyms:	SELV
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)