

Product datasheet for **GA117681**

Human SHISAL2A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SHISAL2A CRISPRa kit - CRISPR gene activation of human shisa like 2A
Symbol:	SHISAL2A
Components:	GA117681G1, SHISAL2A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117681G2, SHISAL2A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117681G3, SHISAL2A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6UWV7</u>
Synonyms:	PRO7171; WWLS2783
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)