

Product datasheet for **GA117692**

Human NUP43 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	NUP43 CRISPRa kit - CRISPR gene activation of human nucleoporin 43
Symbol:	NUP43
Components:	GA117692G1, NUP43 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117692G2, NUP43 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117692G3, NUP43 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8NFH3</u>
Synonyms:	bA350J20.1; FLJ13287; nucleoporin 43kDa; nucleoporin Nup43; p42; p42, FLJ13287, bA350J20.1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)