

Product datasheet for **GA117709**

Human LCE1C activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LCE1C CRISPRa kit - CRISPR gene activation of human late cornified envelope 1C
Symbol:	LCE1C
Components:	GA117709G1, LCE1C gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117709G2, LCE1C gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117709G3, LCE1C gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5T751</u>
Synonyms:	LEP3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)