

Product datasheet for **GA117715**

Human LCE2D activation kit by CRISPRa

Product data:

| | |
|----------------------|--|
| Product Type: | CRISPRa Kits |
| Description: | LCE2D CRISPRa kit - CRISPR gene activation of human late cornified envelope 2D |
| Symbol: | LCE2D |
| Components: | GA117715G1, LCE2D gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117715G2, LCE2D gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117715G3, LCE2D gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID: | <u>Q5TA82</u> |
| Synonyms: | LEP12; SPRL1A |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:

