

Product datasheet for **GA117719**

Human LCE3E activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LCE3E CRISPRa kit - CRISPR gene activation of human late cornified envelope 3E
Symbol:	LCE3E
Components:	GA117719G1, LCE3E gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117719G2, LCE3E gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117719G3, LCE3E gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5T5B0</u>
Synonyms:	LEP17
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)