

Product datasheet for **GA117726**

Human ZNF445 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF445 CRISPRa kit - CRISPR gene activation of human zinc finger protein 445
Symbol:	ZNF445
Components:	GA117726G1, ZNF445 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117726G2, ZNF445 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117726G3, ZNF445 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P59923</u>
Synonyms:	ZKSCAN15; ZNF168; ZSCAN47
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)