

Product datasheet for **GA117729**

Human ANKRD37 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ANKRD37 CRISPRa kit - CRISPR gene activation of human ankyrin repeat domain 37
Symbol:	ANKRD37
Components:	GA117729G1, ANKRD37 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117729G2, ANKRD37 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117729G3, ANKRD37 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q7Z713</u>
Synonyms:	Lrp2bp
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)