

## Product datasheet for GA117742

### Human LPLUNC3 (BPIFB3) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	BPIFB3 CRISPRa kit - CRISPR gene activation of human BPI fold containing family B member 3
Symbol:	BPIFB3
UniProt ID:	<a href="#">P59826</a>
Synonyms:	C20orf185; LPLUNC3; RYA3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)