

Product datasheet for **GA117748**

Human Polyserase1 (TMPRSS9) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMPRSS9 CRISPRa kit - CRISPR gene activation of human transmembrane serine protease 9
Symbol:	TMPRSS9
Components:	GA117748G1, Polyserase1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117748G2, Polyserase1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117748G3, Polyserase1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q7Z410</u>
Synonyms:	FLJ16193
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)