

Product datasheet for **GA117758**

Human **PTCHD3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	PTCHD3 CRISPRa kit - CRISPR gene activation of human patched domain containing 3
Symbol:	PTCHD3
Components:	GA117758G1, PTCHD3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117758G2, PTCHD3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117758G3, PTCHD3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3KNS1</u>
Synonyms:	PTR
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)