


Product datasheet for **GA117768**

Human PTPRQ activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PTPRQ CRISPRa kit - CRISPR gene activation of human protein tyrosine phosphatase receptor type Q
Symbol:	PTPRQ
Components:	GA117768G1, PTPRQ gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117768G2, PTPRQ gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117768G3, PTPRQ gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	DFNB84; DFNB84A; PTPGMC1; R-PTP-Q
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)